



# micomme i series

COPD solution P1 BiPAP



Preheat



Separable  
design



High Technique Perfect Service


## Parameters

|                         |                                                                                     |
|-------------------------|-------------------------------------------------------------------------------------|
| Model                   | PI                                                                                  |
| Mode                    | CPAP, S, T      S/T, VAT, PC                                                        |
| Pressure range          | 4-30cm H <sub>2</sub> O                                                             |
| Pressure accuracy       | ± 0.2 cmH <sub>2</sub> O                                                            |
| Max. operation pressure | 30 cmH <sub>2</sub> O                                                               |
| COMF pressure relief    | 1-3 levels                                                                          |
| Ramp time               | 0 to 45 min (5-minute increments)                                                   |
| Humidification level    | 1-5 levels ( 113 to 185°F/23 to 85°C)                                               |
| Rise time               | 1-6 levels (S,T,S/T)                                                                |
| Data storage capacity   | 8G USB disk                                                                         |
| Weight                  | 1.72 kg                                                                             |
| Mean sound level        | ≤ 30dB                                                                              |
| Application             | OSA, OSAHS, COPD, ARDS and respiratory failure patient in home, clinics or hospital |

## Safe


### AST

Automatic synchronization technology



### VAT

Volume Assured Technology



### APRT


Automatic Pressure Regulation Technology



## Comfortable

### COMF

Pressure relief at the end of exhalation




### One button control



### Preheat

5 minutes from 23°C to 45 °C



## User-friendly

### Rise Time



### Separable design



### USB storage



Micomme Medical Technology Development Co., Ltd.

+86-4000-2000-33      +86-0731-89916333

intl@micomme.com      en.micomme.com

No.8, Lujing Road, High-Tech Zone, Yuelu Changsha, China

